

ColonAiQ® Clinical Evidence Package

CE Marked, China NMPA registered and Singapore HSA-cleared ColonAiQ®

Compiled on August 14, 2025

ColonAiQ® Clinical Evidence Package

Executive Summary

ColonAiQ® represents a breakthrough in colorectal cancer (CRC) detection and monitoring, utilizing a multi-gene methylation approach targeting circulating tumor DNA (ctDNA) in blood samples. This comprehensive evidence package demonstrates the clinical utility of ColonAiQ® across multiple applications: early detection, molecular residual disease (MRD) monitoring, and recurrence prediction.

****Key Performance Metrics:****

- ****Sensitivity****: 86.1% for CRC detection (all stages)
- ****Specificity****: 91.9% for healthy controls
- ****Early Stage Detection****: 78.3% sensitivity for Stage I CRC
- ****MRD Detection****: 78.0% sensitivity with 90.2% specificity at postoperative month 1
- ****Recurrence Prediction****: 17.5x higher risk of relapse for ctDNA-positive patients

1. Technology Overview

Multi-Gene Methylation Platform

ColonAiQ® employs a simplified multiplex quantitative PCR (qPCR) assay targeting six DNA methylation markers:

****Target Genes:****

- ****SEPTIN9**** (2 loci): Well-established CRC methylation marker
- ****BCAT1****: Branched-chain amino acid transaminase 1
- ****IKZF1****: Ikaros family zinc finger 1
- ****VAV3****: Vav guanine nucleotide exchange factor 3
- ****BCAN****: Brevican

Technical Advantages

- 1 **Multi-locus Detection**: Captures fuller spectrum of tumor epigenetic heterogeneity
- 2 **Co-methylation Patterns**: Analyzes adjacent CpG sites from the same DNA molecule
- 3 **Simplified Workflow**: Single-reaction PCR-based assay deployable in clinical settings
- 4 **Cost-Effective**: More accessible than next-generation sequencing approaches
- 5 **Rapid Turnaround**: Results available within days, not weeks

2. Pivotal Clinical Studies

2.1 Gastroenterology Study (2021) - Early Detection

- Study Design**: Multicenter case-control study
- Population**: 507 plasma samples from multiple centers
- Primary Endpoint**: CRC and advanced adenoma detection

Key Findings:

Category	Sensitivity	Specificity	Sample Size
All CRC	86.1%	91.9%	173 CRC, 136 controls
Stage I	78.3%	-	23 patients
Stage II	82.0%	-	50 patients
Stage III	86.0%	-	72 patients
Stage IV	100%	-	16 patients
Advanced Adenoma	42.1%	-	107 patients

Comparative Performance:

- vs. Fecal Immunochemical Test (FIT)**: 88.3% vs 59.7% overall sensitivity
- vs. CEA**: 80-92% vs 17-47% for early-stage CRC
- vs. SEPT9 alone**: Superior performance across all stages

2.2 JAMA Oncology Study (2023) - Molecular Residual Disease

- Study Design**: Multicenter prospective longitudinal cohort study
- Population**: 299 patients with stage I-III CRC

****Study Period****: December 2019 - February 2022

****Follow-up****: Median 21 months (range 8-27 months)

****Primary Findings.****

Postoperative Month 1 Results

- ****ctDNA-positive patients****: 17.5x higher relapse risk (HR 17.5; 95% CI: 8.9-34.4; P<0.001)
- ****Sensitivity****: 78.0% for relapse detection
- ****Specificity****: 90.2%
- ****Negative Predictive Value****: 94.4%

Longitudinal Monitoring

- ****Median Lead Time****: 3.3 months (IQR: 0.5-6.5 months) ahead of radiological detection
- ****Sustained ctDNA-positive****: 68.8x higher recurrence risk (HR 68.8; 95% CI: 18.4-257.7; P<0.001)
- ****Post-adjuvant chemotherapy****: 13.8x higher recurrence risk for ctDNA-positive patients

Clinical Risk Stratification

****Stage III Patients:****

- ctDNA status superior to traditional clinical risk factors
- May guide adjuvant chemotherapy duration decisions
- 6-month vs 3-month chemotherapy benefit observed in high-risk ctDNA-positive patients

2.3 Performance Evaluation Study (Alternative Therapies)

****Study Design****: Performance validation study

****Population****: 88 participants (47 CRC patients, 41 controls)

****Objective****: Comprehensive kit performance assessment

****Performance Metrics.****

- ****Overall Accuracy****: 93.18% (82/88 samples)
- ****CRC Sensitivity****: 89.36% (42/47 patients)
- Stage I: 100% (8/8)
- Stage II: 81.25% (13/16)
- Stage III: 87.5% (14/16)
- Stage IV: 100% (7/7)
- ****Specificity****: 97.56% (40/41 controls)

****Technical Validation:****

- ****Precision****: 100% reproducibility (intra-day and inter-day)
- ****Detection Limit****: 0.5% methylated DNA
- ****Minimum DNA Input****: 10 ng
- ****Interference****: No effect from common clinical substances

2.4 PreC Real-World Study (Annals of Oncology 2023)

****Study Design****: Prospective community screening study

****Population****: 105,285 participants (ages 40-80)

****Study Period****: January 2021 - December 2022

****Setting****: Community-based screening program

****Screening Results:****

- ****ColonAiQ® Positive Rate****: 6.42% (6,759/105,285)
- ****Colonoscopy Compliance****: 48.56% (3,282/6,759) vs 17.25% national average
- ****Positive Predictive Value****:
- CRC: 1.92% (63/3,282)
- Advanced Adenoma: 13.44% (441/3,282)
- Any Adenoma: 36.41%
- ****Early Detection Rate****: 90.28% for early colorectal neoplasms

3. Regulatory Status and Approvals

3.1 China NMPA Approval

- ****Registration Number****: ■■■■20243400902
- ****Approval Date****: May 14, 2024
- ****Valid Until****: May 13, 2029
- ****Classification****: Class C IVD
- ****Intended Use****: Qualitative detection of gene methylation in human plasma samples for CRC auxiliary diagnosis and recurrence risk assessment

3.2 Singapore HSA Registration

- ****Device Registration****: DE0510590

- ****Listing Date****: November 4, 2025
- ****Expiry Date****: October 4, 2026
- ****Classification****: Class C IVD
- ****Quality System****: ISO13485 certified (expires March 9, 2028)

3.3 European CE Marking

- ****Registration****: NL-CA002-2021-63348
- ****Notification Date****: December 6, 2021
- ****Authorized Representative****: SUNGO Europe B.V.
- ****Compliance****: EU Directive 98/79/EG for IVD medical devices

4. Clinical Applications and Utility

4.1 Primary Screening

****Target Population****: Asymptomatic individuals aged 40-80

****Clinical Scenario****: Alternative to invasive colonoscopy for initial screening

****Benefits****:

- Non-invasive blood test
- Higher patient compliance
- Cost-effective first-line screening
- Effective risk stratification

4.2 Molecular Residual Disease Detection

****Target Population****: Post-surgical CRC patients (stages I-III)

****Clinical Scenario****: Early detection of microscopic residual disease

****Optimal Timing****: Postoperative month 1

****Clinical Impact****:

- Guides adjuvant chemotherapy decisions
- Identifies high-risk patients requiring intensive monitoring
- Enables personalized treatment approaches

4.3 Recurrence Monitoring

****Target Population****: CRC survivors under surveillance

****Clinical Scenario****: Serial monitoring during follow-up

****Benefits****:

- Earlier detection than imaging (3.3-month lead time)
- Guides frequency of radiological examinations
- Enables early intervention strategies

4.4 Treatment Response Assessment

****Target Population****: Patients receiving adjuvant chemotherapy

****Clinical Scenario****: Monitoring treatment effectiveness

****Applications****:

- Duration of adjuvant therapy guidance
- Treatment escalation decisions
- Response assessment

5. Proposed Clinical Decision Framework

5.1 ctDNA-Guided Management Protocol

****Positive ctDNA at POM1****

- 1 Recommend adjuvant chemotherapy regardless of stage
- 2 Consider extended duration for persistently positive patients
- 3 Increase surveillance frequency
- 4 Earlier radiological follow-up

****Negative ctDNA at POM1****

- 1 Potential de-escalation of standard adjuvant chemotherapy (research setting)
- 2 Standard surveillance protocol
- 3 Continued serial ctDNA monitoring

****During Surveillance****

- 1 Positive ctDNA → Intensify imaging studies
- 2 Persistent positivity → Consider early intervention
- 3 Serial monitoring every 3 months for up to 2 years

5.2 Integration with Standard Care

****Complementary to CEA Testing:****

- Combined approach improves diagnostic accuracy (AUC: 0.849)
- Addresses CEA limitations (65.9% of patients CEA-negative preoperatively)
- Enhanced sensitivity (83.3%) and specificity (86.5%) when combined

****Colonoscopy Guidance:****

- Positive ColonAiQ® → Prioritize colonoscopy
- Improved compliance rates (48.56% vs 17.25% national average)
- Cost-effective two-step screening approach

6. Comparative Performance Analysis

6.1 vs. Current Screening Methods

Method	Sensitivity (CRC)	Specificity	Notes
ColonAiQ®	86.1%	91.9%	All stages, blood-based
FIT	59.7%	-	Stool-based, affected by bleeding
CEA	17-47%	-	Limited early-stage sensitivity
SEPT9 alone	Lower	-	Single-marker limitation
Colonoscopy	~95%	~95%	Gold standard but invasive

6.2 vs. Other ctDNA Methods

****Advantages over NGS-based approaches:****

- Simplified workflow (hours vs days)
- Lower cost and complexity
- Broader clinical applicability
- Standardized methylation targets
- Reduced technical expertise requirements

7. Health Economic Considerations

7.1 Cost-Effectiveness Factors

****Screening Benefits:****

- Reduced colonoscopy burden for negative tests
- Earlier detection leading to better outcomes
- Improved patient compliance and participation

****MRD Monitoring Value:****

- Prevents unnecessary chemotherapy in low-risk patients
- Enables early intervention in high-risk patients
- Reduces long-term healthcare costs through prevention

****Healthcare System Impact:****

- Streamlined patient pathways
- Resource optimization
- Improved population health outcomes

7.2 Real-World Implementation

****PreC Study Economics:****

- Large-scale feasibility demonstrated (>100,000 participants)
- Sustainable compliance rates
- Effective risk stratification in community settings
- Integration with existing healthcare infrastructure

8. Future Development and Research

8.1 Ongoing Clinical Trials

- Interventional trials incorporating ctDNA-guided management
- Larger prospective validation studies
- Health economic outcome studies

8.2 Potential Expansions

- Additional cancer types using similar methylation approach
- Integration with artificial intelligence for enhanced prediction
- Point-of-care testing development

8.3 Personalized Medicine Integration

- Companion diagnostic development
- Treatment selection algorithms
- Risk-adapted surveillance protocols

9. Conclusions

ColonAiQ® represents a significant advancement in colorectal cancer management, offering clinically validated performance across the entire disease spectrum from screening through surveillance. The comprehensive clinical evidence demonstrates:

- 1 ****Robust Early Detection****: Superior performance compared to existing non-invasive methods
- 2 ****Effective MRD Monitoring****: Clinically actionable results with significant prognostic value
- 3 ****Practical Implementation****: Simplified workflow suitable for routine clinical practice
- 4 ****Regulatory Support****: Multiple international approvals confirming safety and efficacy
- 5 ****Real-World Validation****: Large-scale community implementation demonstrating feasibility

The multi-gene methylation approach addresses key limitations of single-marker assays while maintaining the practical advantages of PCR-based detection. With proven clinical utility across screening, MRD detection, and surveillance applications, ColonAiQ® is positioned to transform colorectal cancer management and improve patient outcomes.

References and Supporting Documents

- 1 Cai G, et al. A Multilocus Blood-Based Assay Targeting Circulating Tumor DNA Methylation Enables Early Detection and Early Relapse Prediction of Colorectal Cancer. *Gastroenterology*. 2021;161(6):2053-2056.
- 1 Mo S, et al. Early Detection of Molecular Residual Disease and Risk Stratification for Stage I to III Colorectal Cancer via Circulating Tumor DNA Methylation. *JAMA Oncol*. 2023;9(6):770-778.

- 1 Liu W, et al. Performance Evaluation of the Multi-gene Methylation Detection Kit for Colorectal Cancer. Alternative Therapies. [E-pub ahead of print].
- 1 Ding Y, et al. Two-year update of the prospective evaluation of ColonAiQ (PreC) study. Annals of Oncology. 2023;34(S2):S426.
- 1 NMPA Medical Device Registration Certificate. Registration No.: ■■■■■20243400902. Approved May 14, 2024.
- 1 HSA Singapore Medical Device Registration. Device Registration No.: DE0510590. Listed November 4, 2025.
- 1 European CE Marking Notification. Registration No.: NL-CA002-2021-63348. December 6, 2021.

This clinical evidence package was compiled on August 14, 2025, based on peer-reviewed publications, regulatory approvals, and clinical study data available at the time of preparation.

Note: This document summarises published and supplied data as of the compilation date. For external communications, always verify the latest official regulatory listings.